


Live Video

Live Video provides a live view of the room, along with the capability to record and bookmark video from an SP or simulation session. Select a camera and then use the camera controls to pan/tilt/zoom the camera view as well as use the customized camera Presets defined in the AV Control application. Bookmarks can also be grouped in accordance with core study areas for the learner such as patient history, CPR, and areas for review.

IMPORTANT!

- Access to Live video is available by permission – consult your administrator.
- If Security Access is enabled and you do not have access to a view based on Room/Organization Unit restrictions, any room you are not allowed to access will be grayed out on the landing page grid when accessing SIMULATIONiQ Enterprise through Internet Explorer.



Home / Video Add New

Explore Search

- Video
 - Live
 - Playback
 - Playlist
 - Pending
 - Bookmarks
 - Video Wall

| <input type="checkbox"/> | Room | Description | Schedule Details |
|--------------------------|------------------------------------|--------------------------------------|------------------|
| <input type="checkbox"/> | ROOM 01- 4220 Cams | Rooms connected with 2, 4220 Cameras | |
| <input type="checkbox"/> | ROOM 02- 5220 Cams | Room connected with 2, 5220 cameras | |
| <input type="checkbox"/> | ROOM 03- 527A Cams | Rooms Connected with 1, 527A Cams | |
| <input type="checkbox"/> | ROOM 04-TVI DVR | Room connected with TVI DVR | |
| <input type="checkbox"/> | ROOM 05- HDSDI DVR | ROOM connected with HDSDI DVR | |
| <input type="checkbox"/> | ROOM 06- NEW DVR | Room Connected with NEW DVR | |

6 record(s) found... 1 1


1. Click **Video > Live Video** to view live video recording as well as record and bookmark video.

Note: Live viewing requires proper viewing privileges designated by your Administrator.

2. Click the desired room to view live video for that room.

Note: If Security Access is enabled and a user does not have access to a view based on Room/Organization Unit restrictions, any room that a user is not allowed to access will be grayed out on the landing page grid (if accessing SIMULATIONiQ Enterprise through Internet Explorer).

Select Camera



Within the web video application, click the link to the predefined camera selection(s) you wish to control for live viewing (as shown above).

Camera Controls



1. Use the following camera preset positions (from **Camera Control**). Use the Pan, Tilt and Zoom In/Out controls for camera placement. A maximum of four presets are allowed for each camera.
2. From the **General Information** area, press **Start Recording** (the **Not Recording** selection appears).
3. Click **Stop Recording** to end the video (the **Recording** selection appears).
4. Click **Refresh** to restore the original screen settings; press to transmit video through a receiver (you can also specify the length of time you wish to display the live view and minimize or maximize the view).

Note: Camera controls are optional, where available.